# LYFOCULTS® PLUS MICROORGANISMS

For quality control procedures.

## SUMMARY AND EXPLANATION

Lyophilization of microorganisms for maintenance and preservation has improved in recent years. Initially, some methods utilized suspensions in fluids such as saline, serum or blood placed in a vacuum with dehydrating agents. Adding gelatin and ascorbic acid was found to be beneficial in the maintenance of the organism during drying.<sup>1</sup> Glucose was found to have a marked preservative effect independent of the drying process employed.<sup>2</sup> The organism pellicle employs a 0.1 ml tryptic soy broth suspension of the appropriate organism with the addition of gelatin, skim milk, ascorbic acid and dextrose to the suspending medium.<sup>3</sup> LyfoCults® Plus uses a modification of this formulation to produce an easy-to-use, cost effective and stable cell preparation.

LyfoCults Plus is a ready-to-use vial containing stabilized, lyophilized, viable microorganisms in an easy-to-use format. These reference organisms are intended for use in laboratory quality control procedures.

## PRINCIPLE

Microbial cultures preserved by lyophilization retain their viability. Ingredients used in preparing the suspensions preserve and protect the microorganisms for use when needed. Following the recommended rehydration procedure, they are immediately ready to use.

## CONTENT

Refer to [www.biomerieux.com](http://www.biomerieux.com) for a complete list of LyfoCults products.

## REAGENTS\*

For *in vitro* diagnostic use.

### LyfoCults Plus

**Active Ingredient:** Microorganism

### Cryoprotectant:

|           |                   |           |
|-----------|-------------------|-----------|
| Albumin   | Ascorbic Acid     | Dextrose  |
| Glycerol  | Mannitol          | Gelatin   |
| Skim Milk | Tryptic Soy Broth | Trehalose |

Rehydration ampoule: 0.6 ml, 0.45% saline

\*Adjusted as required to meet performance standards.

**CAUTION:** LyfoCults Plus contains viable microorganisms and is potentially infectious. All containers and subsequent cultures should be handled according to appropriate CDC/NIH Biosafety Level guidelines and disposed of according to applicable health, safety and environmental regulations.<sup>4</sup>

### Additional materials required

Standard microbiological supplies and equipment such as those commonly found in a microbiology laboratory are not provided.

### Storage instructions


LyfoCults Plus should be stored in its original unopened container at 2-8°C until used or until expiration. When stored properly, organisms will retain viability until the expiration date. **Do not freeze or reuse LyfoCults Plus after initial rehydration. bioMérieux does not warranty the product if reuse occurs or if used after the expiration date.**

### Chemical or physical indications of instability

Each LyfoCults Plus should contain a beige pellicle at the bottom of the container. This product should not be used if there are any signs of deterioration such as premature rehydration, color change, or damage to the container or outer package.

**LyfoCults Plus Rehydration Procedure:**

1. Allow the LyfoCults Plus device to come to room temperature in the pouch. Tear open foil pouch at notch and remove LyfoCults Plus device. (Fig. 1)
2. Remove label from pouch and place on the plate containing the culture media. Remove the protective wrap using the red tear strip. (Fig. 2)
3. Carefully unscrew applicator from vial to vent, and activate ampoule by squeezing button firmly. (Fig. 3)
4. With the applicator at a slight angle to the vial, carefully and slowly squeeze the ampoule dispensing the hydration fluid into the vial. (Fig. 4)
5. Re-insert applicator into vial and gently mix from side-to-side. **Do not invert.** Loosen the ampoule device and hold the vial at a slight angle. Aspirate and expel the suspension a few times with the ampoule device, until the suspension appears homogeneous. The suspension must be used within 30 minutes of rehydration. (Fig. 5)
6. Hold the ampoule device vertically and place 2 drops of resulting suspension on the culture medium. (Fig. 6)
7. Spread the suspension over the medium. For convenience, the tip of the device may be used for isolation streaking or, alternately, a loop may be used. Incubate appropriately (see Table 1). (Fig. 7)
8. Discard used LyfoCults Plus in an appropriate biohazard receptacle. (Fig. 8)

**QUALITY CONTROL**

Each lot number has been tested and approved using standard microbiological techniques and procedures. If the rehydrated culture fails to yield the expected results for the specific reference microorganism strain, the product should not be used for any quality control tests or procedures.

**TABLE 1 - SUGGESTED CULTURE MEDIA AND INCUBATION CONDITIONS**

| Organism                                        | Medium                            | Inoculation | Incubations                                    |
|-------------------------------------------------|-----------------------------------|-------------|------------------------------------------------|
| Nonfastidious Aerobes and Facultative Organisms | TSA with 5% Sheep Blood           | 100 ul      | Aerobic<br>24-48 hours @ 35-37°C               |
| CO <sub>2</sub> Dependent Organisms             | Chocolate Agar                    | 100 ul      | 5-10% CO <sub>2</sub><br>24-72 hours @ 35-37°C |
| Anaerobic Organisms                             | Columbia Agar                     | 100 ul      | Anaerobic<br>48-72 hours @ 35-37°C             |
| <i>Campylobacter</i> spp.                       | Brucella Agar with 5% Sheep Blood | 100 ul      | Micro-aerobic<br>72 hours @ 35-37°C            |
| Fungi: Yeasts and Molds                         | Sabouraud Dextrose Agar           | 100 ul      | Aerobic<br>2-5 days @ 25-30°C                  |
| <i>Legionella</i> spp.                          | Buffered CYE Agar                 | 100 ul      | Aerobic<br>2-5 days @ 35-37°C                  |
| <i>Bordetella</i> spp.                          | TSA with 5% Sheep Blood           | 100 ul      | Aerobic<br>2-4 days @ 35-37°C                  |

**LIMITATIONS OF THE TEST**

LyfoCults Plus contains a high-titered, non-quantitative amount of microorganisms. The rehydration suspension is not intended for use in quantitative analysis.

**REFERENCES**

1. Stamp, Lord. The preservation of bacteria by drying. J Gen Microbiol, 1947, 1:251-265.
2. Annear, Austral. Observations on drying bacteria from the frozen and from the liquid state. J Exp Biol, 1958, 36:211-222.
3. Obara YS, Yamai T, Nikkawa Y, Shimoda, Miyamoto Y. Preservation and transportation of bacteria by a simple gelatin disk method. J Clin Microbiol, 1981, 14:61-66.
4. Richmond JY, McKinney RW. CDC/NIH Manual, Biosafety in Microbiological and Biomedical Laboratories 4th Edition, US Department of Health and Human Services, Washington, D.C., 1999.

**AVAILABILITY**

For technical assistance in the USA, contact bioMérieux Customer Service at 1-800-682-2666. Outside the USA, contact your local bioMérieux Representative.

bioMérieux, the blue logo, and LyfoCults are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

Any other name or trademark is the property of its respective owner.

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

©BIOMÉRIEUX 2010

**INDEX OF SYMBOLS**

| Symbol | Meaning                                             |
|--------|-----------------------------------------------------|
|        | Catalog number                                      |
|        | For <i>in vitro</i> diagnostic use                  |
|        | Manufacturer                                        |
|        | Temperature limitation                              |
|        | Use by                                              |
|        | Batch code                                          |
|        | Consult Instructions for Use                        |
|        | Contains sufficient for <n> tests                   |
|        | Authorized Representative in the European Community |
|        | Biological Risks                                    |

**WARRANTY**

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

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**REFERENCES BIBLIOGRAPHIQUES**

1. Stamp, Lord. The preservation of bacteria by drying. J Gen Microbiol 1947 ; 130 : 1871-82.
2. Annear, Austral. Observations on drying bacteria from the frozen and from the liquid state. J Exp Biol, 1958, 36:211-222.
3. Obara YS, Yamai T, Nikkawa Y, Shimoda, Miyamoto Y. Preservation and transportation of bacteria by a simple gelatin disk method. J Clin Microbiol 1981 ; 5 : 439.
4. Richmond JY, McKinney RW. CDC/NIH Manual, Biosafety in Microbiological and Biomedical Laboratories 4th Edition, US Department of Health and Human Services, Washington, D.C., 1999.

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1. Stamp, Lord. The preservation of bacteria by drying. J Gen Microbiol, 1947, 1:251-265.
2. Annear, Austral. Observations on drying bacteria from the frozen and from the liquid state. J Exp Biol, 1958, 36:211-222.
3. Obara YS. Yamai T, Nikkawa Y, Shimoda, Miyamoto Y. Preservation and transportation of bacteria by a simple gelatin disk method. J Clin Microbiol, 1981, 14:61-66.
4. Richmond JY, McKinney RW. CDC/NIH Manual, Biosafety in Microbiological and Biomedical Laboratories 4th Edition, US Department of Health and Human Services, Washington, D.C., 1999.

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**WARRANTY (Garantía)**

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMÉRIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMÉRIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

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bioMérieux, S.A.