

TRYPTIC SOY MEDIA (SOYBEAN-CASEIN DIGEST MEDIA)

PRODUCTS:

Plated Media:^a

Tryptic Soy Agar	P2500, P3552 (contact plate)
Tryptic Soy Agar with Lecithin and Polysorbate 80	P3500 (contact plate)
Tryptic Soy Agar with Sheep Blood	P2600, P8893 (150mm)
Tryptic Soy Yeast with Rabbit Blood	P2560

Bottled and Tubed Media:^a

Tryptic Soy Agar ^b	B8800, B8805, T8720
Tryptic Soy Broth	Item no. varies by volume
Tryptic Soy Yeast Broth	T3710, T7750

^asee catalog for more ordering options

^bindustrial products

PURPOSE:

Tryptic Soy is a general-purpose base media that is used for the cultivation of fastidious microorganisms. It is a Soybean-Casein Digest media and meets the U.S. Pharmacopeia (USP) standards for use in microbiological tests. Slanted tube media, with or without blood, are used for storage and/or shipment of cultures.

PRINCIPLE:

Tryptic Soy is a highly nutritious media and is commonly used as a base media for the cultivation of microorganisms. The addition of blood to Tryptic Soy Agar (TSA) facilitates the growth of more fastidious microorganisms and allows for the interpretation of hemolytic reactions. Hemolytic reactions can be visualized by holding the agar plate in front of a light source or by observing the area surrounding the colony under low power microscopic magnification (10X). Stabbing a loop into the agar several times at the time of media inoculation⁴ can enhance hemolysis.

HEMOLYTIC REACTIONS:

Beta Hemolysis Clear zone surrounding the colony (Complete lysis of the red blood cells).

Alpha Hemolysis..... Green discoloration (Partial lysis of the red blood cells).

Alpha Prime Hemolysis Green discoloration followed by a clear zone.

Gamma Hemolysis No lysis of the red blood cells surrounding the colony.

Hemolysis may vary with the microorganism and/or with the blood source. Rabbit, horse, and sheep blood are used most frequently. Sheep blood is the generally recommended blood for addition to base agar because beta-hemolytic streptococci show a characteristic clear zone on sheep blood. In addition, *Haemophilus haemolyticus* and *Haemophilus parahaemolyticus* (normal oral flora) appear identical to beta-hemolytic streptococci on horse and rabbit blood and are inhibited on sheep blood.

HEMOLYTIC VARIANCES:

Beta-streptococci:

Groups A,B,C,F,G

Enterococci:**

Non-enterococci:

Listeria species:

Haemophilus species:

H. influenzae

H. parainfluenzae

H. haemolyticus

H. parahaemolyticus

SHEEP

Beta/None

Alpha/None

Alpha/None

Beta/None

No growth

No growth

No growth

No growth

HORSE/RABBIT

Beta/None

Beta/Alpha

Alpha/None

Beta/None

None

None

Beta

Beta

** Some strains of *E. faecalis* may show beta hemolysis.¹

FORMULAS*:

Approximate, per liter of deionized filtered water.

- (1) **Tryptic Soy Agar (TSA):**
 Pancreatic Digest of Casein 15.0 g
 Enzymatic Soy Digest 5.0
 Sodium Chloride 5.0
 Agar 15.0
 Final pH 7.3 ± 0.2 at 25°C
- (2) **Tryptic Soy Agar with Sheep Blood:**
 Same as (1) with various concentrations of Sheep Blood.
- (3) **Tryptic Soy Yeast Agar with 5% Rabbit Blood and 10.0g of Yeast Extract:**
 Same as (1) with 50.0 ml of Rabbit Blood and 10.0g of Yeast Extract.
 Final pH 7.3 ± 0.2 at 25°C
- (4) **Tryptic Soy Broth:**
 Pancreatic Digest of Casein 17.0 g
 Enzymatic Soy Digest 3.0
 Dextrose 2.5
 Sodium Chloride 5.0
 Dipotassium Phosphate 2.5
 Final pH 7.3 ± 0.2 at 25°C
- (5) **Tryptic Soy Yeast Broth:**
 Same as (5) with 10.0 g of Yeast Extract.
 Final pH 7.5 ± 0.2 at 25°C

* Adjustments and/or supplements may be required to meet performance standards.

PRECAUTIONS:*

For *in vitro* diagnostic use. Observe approved biohazard precautions.

Storage: Upon receipt store plated media at 2-8°C and broth media at 2-25°C away from direct light. Media should not be used if there are signs of contamination, deterioration (shrinking, cracking, evaporation, or discoloration), or if the expiration date has passed.

Limitations: Tryptic Soy serves as a nonselective media; biochemical and/or serological testing are necessary for the definitive identification of microorganisms isolated.

Tryptic Soy Media may not support the growth of some fastidious microorganisms and the addition of supplements may be required.

Approximately 2% of group A streptococci may be missed if incubated aerobically unless a provision is made to reduce the oxygen tension. It is recommended that several stabs be made into the blood agar at the time of media inoculation. These group A streptococci produce streptolysin O only, an oxygen-labile hemolysin.³

Incubation in increased CO₂ is less than ideal for determining streptococcal hemolysis because an increased concentration of CO₂ in the presence of oxygen increases the streptococcal production of peroxide which makes the red blood cells resistant to the lysins of beta-streptococci. Preferably, incubation should be in an anaerobic atmosphere,¹ or aerobically with stabs (see above).

The presence of dextrose in Tryptic Soy Broth (TSB) and Tryptic Soy Yeast Broth (TSYB) make these broths unsuitable for maintaining stock cultures; the fermentation of dextrose by the microorganisms may acidify the media and lead to the destruction of microorganisms. These broths are not typically used for direct inoculation of specimens but are used for the growth enhancement of isolated organisms.

PROCEDURE: *

Specimen Collection: Information on specimen collection is found in standard reference material. In general, specimens should be protected from extreme heat and cold and should be delivered to the laboratory without delay. TSB and TSYB media are used to enhance the growth of isolated organisms.

Method of Use: Prior to inoculation, the media should be brought to room temperature. Inoculate according to standard microbiological procedures and on plate media streak inoculum so as to obtain isolated colonies.

Inoculation and incubation procedures can be found in appropriate referenced guidelines such as USP or CLSI.

Incubate under conditions that will permit the growth of microbes. In general, incubate at 35°C for 1-3 days for bacteria and at 22.5°C for 1-5 days for yeasts and fungi. Incubate with adequate moisture in either aerobic, capnophilic, or anaerobic environments, depending on the specific microbes to be cultured.

Interpretation: The presence of growth in broth media is indicated by turbidity. Colonies on agar media will have different morphologies depending on the isolate.

Organism

Streptococcus pyogenes
Streptococcus viridans
Streptococcus pneumoniae
Staphylococcus aureus

Staphylococcus epidermidis

Corynebacteria
Listeria monocytogenes
 Yeasts

Colonial Morphology

Small, beta-hemolytic, transparent to opaque, domed, smooth and entire edge.
 Small, alpha-hemolytic, transparent to opaque, domed, smooth and entire edge.
 Small, alpha-hemolytic, round and mucoid with entire edge.
 Average, ± hemolysis, opaque, circular, smooth, raised, white to golden yellow pigment.
 Average, ± hemolysis, opaque, circular, smooth, raised, usually white to colorless.
 Small, grayish colonies.
 Small, beta-hemolytic, transparent, gray to white.
 Small, white to gray in 48-72 hours.

Materials Required but Not Provided: Standard supplies and equipment commonly found in a microbiological laboratory are not provided.

QUALITY CONTROL:*

Media Used:

TSA Plain

Microorganisms Used (ATCC #):

Escherichia coli (8739)
Streptococcus pyogenes (19615)
Staphylococcus aureus (25923)
Haemophilus influenzae (10211)
Aspergillus niger (16404)[†]
Candida albicans (10231)[†]
Bacillus subtilis (6633)[†]
Pseudomonas aeruginosa (9027)[†]
Staphylococcus aureus (6538)[†]

Expected Results:

Growth
 Growth
 Growth
 Growth with XV disk
 Growth
 Growth
 Growth
 Growth
 Growth

TSA with 5% Sheep Blood

Streptococcus pneumoniae (6305)
Streptococcus pyogenes (19615)
Staphylococcus aureus (25923)
Escherichia coli (25922)

Growth, alpha hemolysis
 Growth, beta hemolysis
 Growth
 Growth

TSA with 5% Rabbit Blood

Streptococcus pneumoniae (6305)
Streptococcus pyogenes (19615)
Staphylococcus aureus (25923)
Escherichia coli (25922)
Haemophilus parahaemolyticus (10014)

Growth, alpha hemolysis
 Growth, beta hemolysis
 Growth
 Growth
 Growth, beta hemolysis

Tryptic Soy Broth

Escherichia coli (25922)
Staphylococcus aureus (25923)

Growth
 Growth

Tryptic Soy Yeast Broth

Candida albicans (10231)
Streptococcus pyogenes (19615)

Growth
 Growth

[†]USP recommended test organism; inoculated with less than 100 colony forming units.

User Quality Control: Check for signs of contamination and deterioration. Tryptic Soy Agar should appear translucent and light tan or straw in color. Tryptic Soy Agar with blood should appear opaque and bright red in color. Tryptic Soy Broth should appear clear and light amber in color.

BIBLIOGRAPHY:

1. Facklam, R. R., *Isolation and Identification of Streptococci*, CDC USHEW, Centers for Disease Control, Atlanta, 1980.
2. Finegold, S. M., and E. J. Baron, *Bailey and Scott's Diagnostic Microbiology*, 7th ed., C. V. Mosby, St. Louis, 1986.
3. Koneman, E. W., et al., *Color Atlas and Textbook of Diagnostic Microbiology*, 3rd ed., J. B. Lippincott, Philadelphia, 1988.
4. Murray, P. R., et al., *Manual of Clinical Microbiology*, 8th ed., American Society for Microbiology, Washington, D. C., 2003.
5. United States Pharmacopeia 30 - NF 25, Chapter 62, Microbial examination of nonsterile products: Tests for specified microorganisms, 2007.

* For more detailed information, consult appropriate references.

PML Microbiologicals, Inc.

Data #780

Copyright 1989 by PML Microbiologicals, Inc.

Revision Date: July 2008